

Shiny Bright不锈钢丝


35.6 cm (14")筒装弓丝

弓丝尺寸	产品代码
0.30mm (.012")	380-238
0.36mm (.014")	380-239
0.41mm (.016")	380-240
0.46mm (.018")	380-241
0.51mm (.020")	380-242
0.41 x 0.41mm (.016" x .016")	380-245
0.41 x 0.56mm (.016" x .022")	380-246
0.43 x 0.64mm (.017" x .025")	380-247
0.46 x 0.56mm (.018" x .022")	380-248
0.46 x 0.64mm (.018" x .025")	380-249
0.48 x 0.48mm (.019" x .019")	380-250
0.48 x 0.64mm (.019" x .025")	380-251
0.51 x 0.56mm (.020" x .022")	380-252
0.53 x 0.64mm (.021" x .025")	380-253
0.56 x 0.71mm (.022" x .028")	380-254

10根/筒

Shiny Bright弓丝具有良好的刚性，同时不失弹性，工艺精湛，尺寸精确，表面光滑，适用于矫正后期的稳定弓形及关闭牙齿间隙。

为方便确认，预成弓丝均标有中线，上颌为三线标记，下颌为单线标记。



弓丝尺寸	平直弓形		方圆形		卵圆形	
	上颌	下颌	上颌	下颌	上颌	下颌
0.30mm (.012")	381-047	381-046	381-049	381-048	383-561	383-560
0.36mm (.014")	381-051	381-052	381-061	381-062	383-563	383-562
0.41mm (.016")	381-053	381-054	381-063	381-064	383-565	383-564
0.46mm (.018")	381-055	381-056	381-065	381-066	383-567	383-566
0.51mm (.020")	381-057	381-058	381-067	381-068	383-569	383-568
0.41 x 0.41mm (.016" x .016")	381-071	381-070	381-081	381-080	383-571	383-570
0.41 x 0.56mm (.016" x .022")	381-073	381-072	381-083	381-082	383-573	383-572
0.43 x 0.64mm (.017" x .025")	381-075	381-074	381-085	381-084	383-575	383-574
0.46 x 0.56mm (.018" x .022")	381-245	381-246	381-091	381-090	383-583	383-582
0.46 x 0.64mm (.018" x .025")	381-077	381-076	381-087	381-086	383-585	383-584
0.48 x 0.64mm (.019" x .025")	381-079	381-078	381-089	381-088	383-577	383-576
0.55 x 0.64mm (.0215" x .025")	381-253	381-254	381-193	381-194	383-579	383-578
0.55 x 0.71mm (.0215" x .028")	381-197	381-198	381-191	381-192	383-581	383-580

10根/包




弓丝尺寸	Bio-Arch I	Bio-Arch II	Bio-Arch III	Bio-Arch IV	Bio-Arch V
0.30mm (.012")	387-001	387-012	387-023	387-034	387-045
0.36mm (.014")	387-002	387-013	387-024	387-035	387-046
0.41mm (.016")	387-003	387-014	387-025	387-036	387-047
0.46mm (.018")	387-004	387-015	387-026	387-037	387-048
0.51mm (.020")	387-005	387-016	387-027	387-038	387-049
0.41 x 0.41mm (.016" x .016")	387-006	387-017	387-028	387-039	
0.41 x 0.56mm (.016" x .022")	387-007	387-018	387-029	387-040	387-051
0.43 x 0.64mm (.017" x .025")	387-008	387-019	387-030	387-041	387-052
0.46 x 0.56mm (.018" x .022")	387-056	387-058	387-060	387-062	387-064
0.46 x 0.64mm (.018" x .025")	387-057	387-059	387-061	387-063	387-065
0.48 x 0.64mm (.019" x .025")	387-009	387-020	387-031	387-042	387-053
0.55 x 0.64mm (.0215" x .025")	387-010	387-021	387-032	387-043	387-054
0.55 x 0.71mm (.0215" x .027")	387-011	387-022	387-033	387-044	387-055

10根/包

三股麻花丝

高度抛光的麻花方丝有较大的工作空间，在治疗早期整平、排齐阶段提供中等大小的力。弓丝切断时不会散开，在弯制时不易折断。为了方便确认，预成弓丝都标有中线，上颌为三线标记，下颌为单线标记。



弓丝尺寸	方圆形	
	上颌	下颌
0.41 x 0.41mm (.016" x .016")	381-133	381-132
0.41 x 0.56mm (.016" x .022")	381-135	381-134
0.46 x 0.64mm (.018" x .025")	381-137	381-136

10根/包



Co-Ax[®]六股麻花丝

这种弓丝可以排齐拥挤和扭转的前牙，中心的一根核心丝保持弓丝的稳定，周围包裹的5根弓丝可增强弹性。该弓丝无需调整即可较长时间提供持续的弹力。表面光滑、缠绕紧密的弓丝可在托槽上自由滑动。预成弓丝上有中线标记，上颌为三线标记，下颌为单线标记。

Co-Ax卷装


0.41mm (.016")	0.46mm (.018")	0.51mm (.020")	0.56mm (.022")
231-356	231-358	231-360	231-362

9 m (30 ft.)/卷

Co-Ax 35.6 cm (14")筒装

0.41mm (.016")	0.46mm (.018")
231-355	231-375

10根/筒





弓丝尺寸	平直弓形		方圆形	
	上颌	下颌	上颌	下颌
0.41mm (.016")	381-127	381-126	381-121	381-120
0.46mm (.018")	381-129	381-128	381-123	381-122

10根/包


八股麻花丝

八股麻花丝具有高弹性和良好的转矩控制，使用范围广。八股麻花丝紧密编织在一起，制成表面光滑的方丝，可以焊接以及弯制成各种圈簧。预成弓丝符合大多数上下颌牙弓的形态。中线标记明显：上颌用三线标记，下颌用单线标记。



弓丝尺寸	平直弓形		方圆形	
	上颌	下颌	上颌	下颌
0.41 x 0.41mm (.016" x .016")	381-203	381-204	381-207	381-208
0.41 x 0.56mm (.016" x .022")	381-151	381-150	381-145	381-144
0.46 x 0.64mm (.018" x .025")	381-153	381-152	381-147	381-146
0.48 x 0.64mm (.019" x .025")	381-201	381-202	381-205	381-206
0.53 x 0.64mm (.021" x .025")	381-155	381-154	381-149	381-148

10根/包



麻花丝是从圆丝到方丝之间很好的过渡弓丝。置入方丝弓槽沟进行转矩控制，所以优良的弹性也用于早期整平。简单的病例仅用麻花丝就可以达到矫治效果。